



C7C CLAW WITH MOTOR:	PART # 264365
C7C CLAW WITHOUT MOTOR:	PART # 264365-LM
7.5HP 3520 RPM MOTOR:	PART # 244266

1	10429	ORIGINAL RELEASE	10/28/2021	NHR
REV	ECD NO.	DESCRIPTION OF CHANGE	DATE	BY
<small>DO NOT SCALE DRAWINGS UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES AND TOLERANCES ARE:</small> XX ± 0.1 XXX ± 0.01 XXXX ± 0.005 .XXX ± .1		DRAWN BY NHR 10/28/2021 CHECKED BY DATE APPROVED BY DATE		
<small>THIRD ANGLE PROJECTION</small> 		TITLE DESCRIPTION VACUUM PUMP, CLAW (C7C) WITH 7.5HP 213TC NEMA MOTOR	COPYRIGHT & COPYRIGHTED BY OHIO MEDICAL LLC	
MODEL N/A	FINISH N/A	Dwg No. 264365	REV 1	SHEET 1 OF 1
FILE NAME 264365.DWG		SCALE: 1:4	CONFIDENTIAL: THIS DRAWING AND THE DESIGN REPRESENTED THEREIN ARE THE PROPERTY OF OHIO MEDICAL LLC AND HAS BEEN ISSUED WITH THE UNDERSTANDING THAT NEITHER DRAWING NOR THE INFORMATION THEREIN IS TO BE USED OR DISCLOSED, EITHER IN WHOLE OR PART, EXCEPT AS SPECIFICALLY AUTHORIZED BY OHIO MEDICAL LLC AND SHALL BE RETURNED UPON REQUEST.	