





ITEM NO.	QTY.	PART NUMBER	DESCRIPTION	REF NOTE
1	2	233907		
2	4	242005	1" SOURCE ISOLATION VALVE	
3	12	A-12105-6	1" NPT CLOSE NIPPLE	
4	2	A-12108-6	1" Threaded Tee	
5	4	A-12109-6	1 NPT - ELBOW	
6	4	A-12137-6	UNION, 1"	
7	4	A-14351-3	1" NPT SCH 40 GALVANIZED PIPE (8" LONG)	
8	1	BASE-MED-FB-3	WELDED BASE, MEDICAL FILTER BANK 3" PIPE	



REV	ECO NO.	DESCRIPTION OF CHANGE	DATE	BY
01	N/A	INITIAL RELEASE	02/05/16	F.L.

DO NOT SCALE DRAWING	DRAWN BY F. LOPEZ	DATE 02/05/16	
UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES AND TOLERANCES ARE: .XX ± 0.050 .XXX ± 0.020 XXXX ± 0.005	CHECKED BY	DATE	
THIRD ANGLE PROJECTION	APPROVED BY	DATE	COPYRIGHT & COPYRIGHT © BY OHIO MEDICAL LLC
MATERIAL N/A	ITEM DESCRIPTION MEDICAL VACUUM FILTER BANK, 1" NPT CONNECTIONS		
FINISH N/A	DWG SIZE B	DRAWING NUMBER BOM-MED-FB-1	REV 01
FILE NAME BOM-MED-FB-1.SLDDRW	SCALE: 1:1		SHEET 1 OF 1

CONFIDENTIAL: THIS DRAWING AND THE DESIGN REPRESENTED THEREIN ARE THE PROPERTY OF OHIO MEDICAL AND HAS BEEN ISSUED WITH THE UNDERSTANDING THAT NEITHER DRAWING NOR THE INFORMATION THEREIN IS TO BE USED OR REPRODUCED, EITHER IN WHOLE OR PART, EXCEPT AS SPECIFICALLY AUTHORIZED BY OHIO MEDICAL AND SHALL BE RETURNED UPON REQUEST.

BOM-MED-FB-1 01