




ITEM NO.	QTY.	PART NUMBER	DESCRIPTION	REF NOTE
1	2	233909	FILTER, INLET MEDICAL VACUUM, 1-1/2" FPT	
2	4	242007	1-1/2" NPT BALL VALVE	
3	4	A-12102-1	1.5 X 3" NPT NIPPLE	
4	4	A-12102-33	1.5 X 2" NPT NIPPLE	
5	4	A-12105-8	1-1/2" NPT CLOSE NIPPLE	
6	2	A-12108-8	1-1/2" Threaded Tee	
7	4	A-12109-8	1-1/2 NPT - ELBOW	
8	4	A-12137-8	UNION, 1-1/2"	
9	4	A-14351-5	1-1/2" NPT SCH 40 GALVANIZED PIPE (6" LONG)	
10	1	BASE-MED FB-3	WELDED BASE, MEDICAL FILTER BANK 3" PIPE	

1-1/2" NPT VACUUM OUTLET CONNECTION

1-1/2" NPT VACUUM INLET CONNECTION



10

01	N/A	INITIAL RELEASE	01/27/17	F.L.
REV	ECO NO.	DESCRIPTION OF CHANGE	DATE	BY
DO NOT SCALE DRAWING		DRAWN BY	DATE	
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES AND TOLERANCES ARE:		F. LOPEZ	01/27/17	
.XX ± 0.1		CHECKED BY	DATE	
.XXX ± 0.01		APPROVED BY	DATE	
X.XXX ± 0.005				
XXX ± 1				
THIRD ANGLE PERSPECTION				
MATERIAL:		ITEM DESCRIPTION:		
N/A		MEDICAL VACUUM FILTER BANK, 1-1/2" NPT CONNECTIONS		
FINISH:		DWG SIZE	DRAWING NUMBER	REV
N/A		B	BOM-MED FB-15	01
FILE NAME:		SCALE:	SHEET	OF
BOM-MED FB-15_SLDRAW		1:1	1	1

CONFIDENTIAL: THIS DRAWING AND THE DESIGN REPRESENTED THEREIN ARE THE PROPERTY OF OHIO MEDICAL AND HAS BEEN ISSUED WITH THE UNDERSTANDING THAT NEITHER DRAWING NOR THE INFORMATION THEREIN IS TO BE USED OR REPRODUCED, EITHER IN WHOLE OR PART, EXCEPT AS SPECIFICALLY AUTHORIZED BY OHIO MEDICAL AND SHALL BE RETURNED UPON REQUEST.