



ITEM NO.	Default/ QTY.	PART NUMBER	DESCRIPTION	REF NOTE
1	2	233913	MEDICAL VACUUM FILTER, 0.2 MICRONS 3" FPT	
2	4	A-12102-34	3 X 4" NPT NIPPLE	
3	4	A-12105-11	3" NPT CLOSE NIPPLE	
4	2	A-12108-11	3" Threaded Tee	
5	4	A-12109-11	3 NPT - ELBOW	
6	4	A-14351-8	3" NPT SCH 40 GALVANIZED PIPE (6" LONG)	
7	8	B-14404-7	FLANGE, 3"	
8	1	BASE-MED FB-3	WELDED BASE, MEDICAL FILTER BANK 3" PIPE	
9	4	BUT-3	BUTTERFLY VALVE 3" NPT	



01	N/A	INITIAL RELEASE	01/14/20	F.L.
REV	ECO NO.	DESCRIPTION OF CHANGE	DATE	BY
DO NOT SCALE DRAWING		DRAWN BY	DATE	
UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES AND TOLERANCES ARE:		F. LOPEZ	01/14/20	
.XX ± 0.050		CHECKED BY	DATE	
.XXX ± 0.020		APPROVED BY	DATE	
XXXX ± 0.005				
THIRD ANGLE PROJECTION		COPYRIGHT & COPYRIGHT D BY OHIO MEDICAL LLC		
MATERIAL	N/A	ITEM DESCRIPTION		
FINISH	N/A	MEDICAL VACUUM FILTER BANK, 3" NPT CONNECTIONS		
FILE NAME	BOM-MED FB-3.SLDDRW	DWG SIZE	DRAWING NUMBER	REV
		B	BOM-MED FB-3	01
			SCALE: 1:1	SHEET 1 OF 1

BOM-MED FB-3 REV 01