




4" NPT VACUUM OUTLET



4" NPT VACUUM INLET

ITEM NO.	QTY.	PART NUMBER	DESCRIPTION	REF NOTE
1	2	233914	MEDICAL VACUUM FILTER, 0.2 MICRON 4" FPT	
2	4	A-12102-41	4 X 5" NPT NIPPLE	
3	4	A-12103-86	4 X 6" NPT NIPPLE	
4	4	A-12105-13	4" NPT CLOSE NIPPLE	
5	2	A-12108-13	4" Threaded Tee	
6	4	A-12109-13	4 NPT - ELBOW	
7	8	B-14404-9	FLANGE, 4"	
8	1	BASE-MED FB-4	WELDED BASE, MEDICAL FILTER BANK 4" PIPE	
9	4	BUT-4	BUTTERFLY VALVE, 4"	



01	N/A	INITIAL RELEASE	04/01/13	F.L.
REV.	ECO NO.	DESCRIPTION OF CHANGE	DATE	BY
DO NOT SCALE DRAWING UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES AND TOLERANCES ARE: X.X ± 0.1 X.XX ± 0.01 X.XXX ± 0.005 X.X" ± 1		DRAWN BY F. LOPEZ	DATE 04/01/13	
		CHECKED BY	DATE	
THIRD ANGLE PROJECTION		APPROVED BY	DATE	COPYRIGHT © AND PATENTED BY OHIO MEDICAL CORPORATION
MATERIAL N/A		ITEM DESCRIPTION MEDICAL VACUUM FILTER BANK, 4" NPT CONNECTIONS		
FINISH N/A		DWG SIZE B	DRAWING NUMBER BOM-MED FB-4	REV. 01
FILE NAME BOM-MED FB-4.slddrw		SCALE: 1:16		SHEET 1 OF 1

REV: 01
BOM-MED FB-4