

Replacement Electrical Control Panels for Air Systems

Ohio Medical electrical control panels are ideal for air system upgrades in hospitals and labs. Equipped with high definition color touch screens, these panels display multiple conditions, including compressor running status, elapsed hours, pressure level, and audio/visual alarms.

Standard features include lag alarm, automatic compressor alternation, adjustable minimum run timer, adjustable pressure settings, adjustable screen brightness, alarm history registry, and data logging.

These panels also monitor dew point and CO levels when fitted with dew point and CO transducers. Or if the facility prefers to use existing dew point and CO monitors, the panel will function without those inputs. The panel can also be provided with 115-volt outputs and switches for dryers.

Connections to building management systems are also included. Alarms are capable of interface via Ethernet, Modbus, etc. The alarms are programmed for both English and Spanish languages.



Use our Easy Configurator at ohiomedicalparts.com/panels to find the panel for your system.