

ULPA Vacuum Filters

Ohio Medical ULPA Vacuum Filters are designed for medical vacuum applications and feature Ultra-Low Particulate Air (ULPA) media rated for 99.999% removal efficiency at 0.1 micron. Each assembly includes a filter housing, Pyrex collection bottle, valve, tubing, and biohazard label. Available with 1-inch through 3-inch NPSC connections, the filters feature replaceable ULPA elements and additional ports for drainage and servicing.



Representative 2-inch assembly shown.

Filter assemblies

Ohio Medical part number	Inlet and outlet connection
233960	1-inch NPSC
233961	1-1/4-inch NPSC
233962	1-1/2-inch NPSC
233963	2-inch NPSC
233964	2-1/2-inch NPSC
233965	3-inch NPSC

Replacement parts

Ohio Medical part number	Description	Fits filter sizes
233970	ULPA replacement element	1, 1-1/4 and 1-1/2 inches
233971	ULPA replacement element	2 and 2-1/2 inches
233972	ULPA replacement element	3 inches
233973	Pyrex replacement bottle with cap	1 through 3 inches